

## IMMEDIATE HEMOGLOBIN RESULTS FROM A HAND-HELD ANALYZER



Results in less than 2  
seconds

Reagent-free  
microcuvettes

40 days or 10,000 tests  
continuous use

No start-up process



### FAST AND ACCURATE

- Results in less than 2 seconds
- Measuring range: 0.3 – 25.5 g/dL
- Imprecision: CV < 1%
- Factory calibrated against HiCN reference method according to ICSH
- Photometric method with broad spectrum and compensation for turbidity and scattering
- No re-calibration required

### EASY TO USE

- 'Always on' technology. The analyzer is always ready for use
- Automatic self-check between every measurement
- Long microcuvette shelf life: 2.5 years from date of production
- Microcuvettes unaffected by moisture or temperature and can be used until expiration date after opening the bag

### PRACTICAL AND PORTABLE

- Sample volume: <10 µL
- Capillary or venous blood (EDTA and Heparin anticoagulant can be used)
- Minimal maintenance
- Operating ambient temperature: 10° – 35° C
- Storage temperature: 0° – 50° C (short term -30° C – 70° C)
- Integrated rechargeable lithium-ion battery (40 days / 10,000 tests of continuous use)
- Compact: 150 (d) x 90 (w) x 40mm (h)
- Lightweight: 185 g
- Connect to PC via USB 2.0 or Bluetooth (optional)

**Fast, easy to use and reliable** in all climate conditions



**1**

Collect blood sample



**2**

Ensure analyzer is ready for use



**3**

Results appear in 1 - 2 seconds

**Distributed by**

**Point of Care Diagnostics Pty Ltd**

Phone: 1800 640 075

Email: [admin@pocd.com.au](mailto:admin@pocd.com.au)

Web: [www.pocd.com.au](http://www.pocd.com.au)

#### **Manufacturer**

**DiaSpect Medical GmbH**

Von-Cancrin-Str. 3

63877 Sailauf

Germany

**Tel.** +49 6093 996689-60

#### **Sales**

**EKF Diagnostics Holdings plc**

Avon House, 19 Stanwell Road,

Penarth, Cardiff, UK, CF64 2EZ

**Tel.** +44 (0)29 20 710570

**[www.ekfdiagnostics.com](http://www.ekfdiagnostics.com)**

Approved for the  
European Market

**EKF**  
**DIAGNOSTICS**

Disclaimer: Diaspect Medical GmbH, an EKF Diagnostics Holdings plc company, reserves the right to change the product specifications without prior notice.

Revision EU 0.0-06/2014