



# Afinion™ HbA1c

Test for quantitative determination of glycated haemoglobin (HbA1c) in human whole blood, used to monitor metabolic control in patients with diabetes. May be used for diagnosis of diabetes and in identifying patients who may be at risk for developing diabetes, when in compliance with national regulations\*.

- **simple** • **fast** • **reliable**
- 3 minutes test time
- Excellent precision
- No interference from Hb variants
- CLIA waived



\*May only be used for monitoring in USA.

## The 3-minute Afinion™ test fits perfectly into the workflow

“Immediate feedback motivates patients to make the necessary treatment changes.”

*Cagliero E et. al. Diabetes Care 1999;22:1785-1789.*

“Use of point-of-care testing for HbA1c allows for timely decisions on therapy changes, when needed.”

*ADA, 2011. Diabetes Care 2011;34(supl 1):s11-s61.*

## 3-step procedure



## Features

- 1.5 µL sample volume
- Configurable display of HbA1c units
- 4-15% HbA1c, 20-140 mmol/mol HbA1c
- 3 months room temperature storage
- Ready to use Controls in 2 levels
- NGSP and IFCC certified

## Other available tests

Afinion™ ACR, Afinion™ CRP and Afinion™ Lipid Panel



### Alere

Australia  
1800 622 642 Free call in AUS  
+61 7 3363 7100 Office  
+61 7 3363 7199 Fax  
au.enquiries@alere.com

[alere.com.au](http://alere.com.au)

### Alere Ltd

New Zealand  
0800 466 951 Free call in NZ  
0800 466 951 Office  
0800 466 958 Fax  
nz.enquiries@alere.com

[alere.co.nz](http://alere.co.nz)

Distributed by: