



Alere™

Afinion™ Lipid Panel

NEW

Test for the quantitative determination of Total Cholesterol (Chol), High-Density Lipoprotein cholesterol (HDL), Low-Density Lipoprotein cholesterol (LDL), Triglycerides (Trig), non-HDL and Chol/HDL ratio in whole blood, serum and plasma, used in the diagnosis and treatment of lipid disorders.

- **simple** • **fast** • **reliable**
- Specially designed sampling device ensures easy and safe sample collection
- Barcode includes calibration and lot specific data - no user calibration necessary
- Fully enclosed testing environment increases user safety



The importance of measuring a lipid panel

It is recommended by the NCEP (National Cholesterol Education Program) that a Lipid Panel/Cholesterol Test should be measured:

- Every 5 years in healthy adults.
- More regularly if you have other risk factors
- In children and youths who are at an increased risk
- To assess the effectiveness of lipid-lowering medication

American Diabetes Association

“In most adult patients, measure fasting lipid profile at least annually. In adults with low-risk lipid values (...) lipid assessments may be tested every 2 years.”

Reference: ADA. Standards of Medical Care in Diabetes - 2011. Diabetes Care 2011;34(supl 1):s11-s61.

3-step procedure



Features

- 15 µL sample volume
- 8 minute test time
- 6 results in one test
- Configurable display of measuring units
- Ready to use control material

Other available tests

Afinion™ HbA1c, Afinion™ CRP and Afinion™ ACR



Alere

Australia
1800 622 642 Free call in AUS
+61 7 3363 7100 Office
+61 7 3363 7199 Fax
au.enquiries@alere.com

alere.com.au

Alere Ltd

New Zealand
0800 466 951 Free call in NZ
0800 466 951 Office
0800 466 958 Fax
nz.enquiries@alere.com

alere.co.nz

Distributed by: