

**Alere**TM

AfinionTM CRP

AfinionTM CRP is a rapid test for quantitative determination of C-reactive protein (CRP) in human blood. CRP measurement is valuable in diagnosing and monitoring of infections and non-infectious inflammatory diseases.

- **simple • fast • reliable**
- 1.5 μ L sample volume
- 4 minutes test time
- Excellent precision
- No user calibration necessary – barcode includes calibration and lot-specific data



Why measure C-reactive protein?

A significantly increased CRP result may indicate the need for immediate antibiotic treatment. The CRP level decreases rapidly in response to effective therapy.

A normal or moderately increased CRP may support a diagnosis of viral or self-limiting bacterial infection, which provides valuable information to prevent unnecessary antibiotic prescriptions.

3-step procedure



Facts

- Accurate whole blood measurement with automatic hematocrit correction
- Excellent correlation with laboratory methods
- Excellent precision in the total measuring range of 5-200 mg/L
- 4 weeks room temperature storage
- Ready to use Controls in 2 levels

Other available tests

Afinion™ HbA1c, Afinion™ ACR, and Afinion™ Lipid Panel



Australia
New Zealand

au.enquiries@alere.com
nz.enquiries@alere.com

1800 622 642 Free call in AUS
0800 466 951 Free call in NZ

alere.com.au
alere.co.nz