

UNIFINE PENTIPS COMPATIBILITY

Unifine Pentips®
NEW
4mm 

0.23mm (32G) x 4mm

Pack Size	100
OM code	AN3540
PIP code	379-6133

Unifine Pentips®
5mm 

0.25mm (31G) x 5mm

Pack Size	100
OM code	AN3550
PIP code	368-0881

Unifine Pentips®
6mm 

0.25mm (31G) x 6mm

Pack Size	100
OM code	AN3590
PIP code	281-5793

Unifine Pentips®
8mm 

0.25mm (31G) x 8mm

Pack Size	100
OM code	AN3530
PIP code	234-8415

Unifine Pentips®
12mm 

0.33mm (29G) x 12mm

Pack Size	100
OM code	AN3529
PIP code	234-8423

Fits **ALL** leading brands of injection pens*



Owen Mumford: Autopen® Classic



Owen Mumford: Autopen® 24

**Available from Owen Mumford
or your usual supplier.
To order call 0800 032 0439.**



Novo Nordisk: Victoza®



Novo Nordisk: NovoPen® 3



Novo Nordisk: NovoRapid® FlexPen®



Novo Nordisk: Levemir® FlexPen®




Novo Nordisk: NovoPen® 4



Novo Nordisk: NovoPen Echo®



Novo Nordisk: NovoMix® FlexPen® 3



Novo Nordisk: NovoPen® Junior



Novo Nordisk: NovoPen® FlexTouch



Lilly: Byetta®



Lilly: Humalog® KwikPen™



Lilly: Humapen® Memoir



Lilly: Humalog® Pen



Lilly: Humapen® Luxura



Lilly: Humulin® Pen



Sanofi: Apidra® SoloStar®



Sanofi: Lantus® SoloStar®



Sanofi: OptiPen® Pro 1 & 2



Sanofi: ClikSTAR®



Sanofi: OptiClik®



Sanofi: OptiSet®

* Compatibility data correct as of 01/10/13.

Unifine® Pentips® and Autopen® are registered trademarks of Owen Mumford Ltd. All rights reserved. All other trademarks are the property of their respective owners.

UNILET LANCET COMPATIBILITY

Unilet® Eco 28G



Gauge: 28G (0.375mm) Tariff: Type A

For maximum comfort

Pack Size	100	200
OM code	AT0935	AT0940
PIP code	357-8903	357-8911

Unilet® ComforTouch™



Gauge: 28G (0.375mm) Tariff: Type A + B

For maximum comfort

Pack Size	100	200
OM code	AT0465	AT0460
PIP code	224-2568	224-2576

Unilet® GP Superlite



Gauge: 23G (0.66mm) Tariff: Type A

Optimised for bloodflow
Ideal for older meters or thicker skin

Pack Size	100	200
OM code	AT0455	AT0450
PIP code	045-2144	045-2169

Compatible with most capillary blood lancing devices*

*Except Roche Multiclix and Softlix Plus



Owen Mumford

Device: Autolet® Impression



Abbott

Device: FreeStyle

Meter: FreeStyle Lite;
FreeStyle Freedom Lite



Bayer

Device: Ascensia MICROLET®

Meter: Ascensia CONTOUR®;
Ascensia BREEZE®



Bayer

Device: Ascensia MICROLET® Vaculance



BBI Healthcare

Device: SensoCard

Meter: SensoCard Plus Talking Meter



BD

Device: Optimus™



Home Diagnostics

Device: Gentle Draw®

Meter: Prestige QX Smart System®;
TRUETrack® Smart System



Lifescan

Device: OneTouch® UltraSoft®

Meter: OneTouch® Ultra®; OneTouch® Ultra®2;
OneTouch® UltraSmart®



Lifescan

Device: Mini OneTouch®

Meter: OneTouch® UltraEasy®



Menarini

Device: Glucoject Dual

Meter: GlucoMen® Visio; GlucoMen® PC

CUSTOMER SERVICES TEAM:
0800 032 0439

Unilet® and Autolet® Impression are registered trademarks of Owen Mumford Ltd. All rights reserved. All other trademarks are the property of their respective owners.